Non-Sterile Atrial Lift System (Rigid)

Catalog Numbers: 400-525, 401-501, 401-502, 401-503, 401-504, 401-505,
401-506, 401-507, 401-508, & 401-509

PRODUCT DESCRIPTION

The Atrial Lift System is reusable system comprised of a Rigid Table Mount Instrument Holder, two Blade Handles (140mm and 200mm lengths) and six Atrial Lift Blades (width and length dimensions of 35 x 40mm, 35 x 50mm, 35 x 60mm, 40 x 50 mm, 40 x 60mm and 40 x 70mm). Blade introduction forceps are included to introduce the atrial lift blade into the surgical incision. The Atrial Lift System is used to lift the atrial wall to provide surgical access to heart structures for surgical procedures using thoracotomy access.

CLEANING

Disassemble and scrub the instruments and parts thoroughly using a soft brush and mild detergent (i.e. Manu-Klenz, manufactured by *Steris Corporation*). Remove all traces of blood and debris. Make sure all moving parts are cleaned thoroughly to prevent debris from interfering with movement. It is recommended that the instruments and parts be ultrasonically cleaned.

It is extremely important that moveable parts be properly lubricated to keep these parts functional. It is recommended that all components be immersed in a water-soluble lubricant (i.e. Metri-Clean, manufactured by *Metrex*).

STERILIZATION

The instruments may be steam or gas sterilized before use and after each cleaning. Refer to the sterilizer manufacturer's instructions for correct time, temperature and pressure settings.

USE OF DEVICE

1. Make appropriate surgical incision to gain access for atrial incision.
2. Once atrial incision is made for surgical repair, determine the size of the atrial lift blade required and place lifter blade in jaws of the introducer forceps. See Figure 1.
3. Determine appropriate location for entry of blade handle shaft through chest wall.
4. Make small puncture incision in chest wall.
5. Mount table mount instrument holder on operating room table (see enclosed table mount instrument holder instructions).
6. Place shaft of blade handle with threads through securing mechanism of Table mount instrument holder and pass shaft of handle through incision in chest wall (see enclosed table mount instrument holder instructions).
7. Using introducer forceps to position lifter blade inside chest cavity so that the male threads of the handle shaft can be threaded into the female threaded receptacle in the lift blade. Thread the shaft into the blade by turning clockwise.
8. Position the lift blade in the atrial incision; once the proper position is achieved tighten the adjustment knobs on the table mount instrument holder (see enclosed table mount instrument holder instructions).
9. To remove lifter, loosen the adjustment knobs on the table mount instrument holder. Grasp lift blade with introducer forceps and remove handle by rotating handle counter-clockwise until shaft is disengaged from blade. Remove blade from chest cavity and withdraw handle shaft from incision (See table mount instrument holder instructions).

WARNINGS

Carefully choose the location of the puncture incision for introduction of the handle shaft and introduce the shaft carefully avoiding vital physiologic structures.

Use care to avoid tissue damage when introducing and lifting with atrial lift blade.

WARRANTY

AtriCure warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond AtriCure control directly affect the instrument and the result obtained from its use. AtriCure's obligation under this warranty is limited to the repair or replacement of this instrument and AtriCure shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. AtriCure neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AtriCure assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.

SERVICE OR REPAIR /TECHNICAL INFORMATION

Contact information for service or repair or to request technical information:

Manufactured by:
AtriCure Incorporated
7555 Innovation Way
Mason, Ohio 45040 USA
Customer Service:
1-866-349-2342 (toll free)
1-513-755-4100 (phone)

European Representative:
Köntges SPRL
Avenue Hellevelt 35
B-1180 Brussels Belgium
Tel: +32 (0) 2 375 51 63
FAX: +32 (0) 2 375 89 06
e-mail: herbert.kontges@skynet.be

Graphic Symbols for Device Labeling


Catalogue Number 	Batch Code 	Quantity 	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner. 
Manufacturer/Company Address 	Authorized Representative in the European Union 	Instrument supplied Non-Sterile 	CE mark, Product conforms with the essential requirements in the European Medical Devices Directive 93/42/EEC. 
Does not contain Natural Rubber Latex 	Date of Manufacture 	Consult Instructions For Use 	Contains no di (2-ethylhexyl) phthalate (DEHP) released from polyvinyl chloride (PVC) 

Figure 1
Atrial Lift Blade and Blade Introduction Forceps

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