AtriCure

AtriClip® Gillinov-Cosgrove® Selection Guide (CGG100)


Instructions for Use

DESCRIPTION

The AtriClip Gillinov-Cosgrove Selection Guide (Guide) is a single patient use surgical instrument designed to assist in the selection of the appropriate AtriClip LAA Exclusion System with preloaded Gillinov-Cosgrove LAA Clip (Clip). The Guide is an accessory to the Clip. The Guide is malleable and may be placed directly adjacent to the Left Atrial Appendage (LAA) for sizing.

There are markings on the guide which facilitate the evaluation of the structure and selection of the appropriate Clip size. The Guide does not contain natural rubber latex components.

AtriClip Gillinov-Cosgrove Selection Guide Illustration

1. Handle
2. Base
3. Indication Marks

⚠ BEFORE USING PRODUCT READ THE FOLLOWING INFORMATION THOROUGHLY

IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

INDICATION FOR USE

The AtriClip Gillinov-Cosgrove Selection Guide (Guide) is intended to be used as an accessory to the AtriClip LAA Exclusion System with preloaded Gillinov-Cosgrove Clip (Clip).

CONTRAINDICATIONS

None known

⚠ WARNINGS

1. Do not use the Guide if the sterility barrier is breached. To avoid the risk of patient infection, inspect the product packaging prior to opening to ensure that the sterility barrier is not breached.
2. Do not apply excessive force when using the Guide. Using excessive force may cause tissue damage.


⚠ PRECAUTIONS

1. Carefully read all instructions for the AtriClip Gillinov-Cosgrove Selection Guide prior to using the device. Failure to properly follow instructions may lead to injury and/or may result in improper use of the device.
2. The AtriClip Gillinov-Cosgrove Selection Guide is to be used only to assist with selecting the appropriate AtriClip LAA Exclusion System with preloaded Gillinov-Cosgrove Clip.
3. Use of the Guide should be limited to properly trained and qualified medical personnel.
4. To avoid damage to the device, do not drop or toss the Guide. If the Guide is dropped, do not use. Replace with a new Guide.
5. To avoid damage to the Guide, use caution when handling and interfacing with other devices.
6. Do not bend the guide in the area of the indication marks.
7. Do not re-sterilize or reuse the Guide. Single Patient use only.

INSTRUCTIONS FOR USE

1. Using sterile technique, remove the Guide from its packaging.
2. Align the Guide adjacent with the LAA.

NOTE: Evaluation should occur at the base of the appendage.

3. Evaluate the LAA and select the appropriate Clip size.

Clip Size	LAA Size Range
35 mm	29 - 35 mm
40 mm	34 - 40 mm
45 mm	39 - 45 mm
50 mm	44 - 50 mm

4. Discard the Guide after use. Follow the local governing ordinances and recycling plans regarding disposal or recycling of device components.

HOW SUPPLIED

The AtriClip Gillinov-Cosgrove Selection Guide is supplied as a STERILE instrument and is for single patient use only.

Sterility is guaranteed unless the package is opened or damaged. Do not re-sterilize. Do not use the Gillinov-Cosgrove Selection Guide after the expiration date indicated.

RETURN OF USED PRODUCT

If for any reason this product must be returned to AtriCure, Inc., a return goods authorization (RGA) number is required from AtriCure, Inc., prior to shipping.

If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from AtriCure, Inc.

⚠ CAUTION: It is the responsibility of the health care institution to adequately prepare and identify the product for shipment.

DISCLAIMER STATEMENTS

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property

EXPLANATION OF SYMBOLS ON PACKAGE LABELING

Refer to the outer package label to see which symbols apply to this product.

	Non-Pyrogenic		Use-by Date
	Sterilized by Irradiation		Caution
	Single Use Only		Consult Instructions for Use
	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner		Manufacturer
	Lot Number		Do not Use if package is damaged
	Temperature Limit		Humidity Limit
	Do Not Resterilize		Authorized Representative in the Brazilian Community

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