
INSTRUCTIONS FOR USE

SYSTEM AND PRODUCT DESCRIPTION

The Isolator® Transpolar Pen, henceforth referred to as the Isolator Pen or the Pen, are used in conjunction with an AtriCure RF Generator (ASU2/ASU3 and ASB3 or MAG) and Footswitch. The Isolator Pen is a sterile, single use electrosurgery device intended to ablate cardiac tissue when connected to a compatible AtriCure radiofrequency Generator. When the Pen is connected to the Generator, the Generator provides bipolar radiofrequency (RF) energy flowing between both electrodes of the Pen. The Operator controls the application of this RF energy by pressing and holding the Footswitch or button. The Isolator Pen may be used for temporary cardiac pacing, sensing, recording, and stimulation during the evaluation of cardiac arrhythmias when connected to a temporary external cardiac pacemaker or recording device.

INDICATION FOR USE

The Isolator Pen is a sterile, single use electrosurgery device intended to ablate cardiac tissue during cardiac surgery using radiofrequency (RF) energy when connected directly to the Generator in Ablation mode.

The Isolator Pen may be used for temporary cardiac pacing, sensing, recording, and stimulation during the evaluation of cardiac arrhythmias during surgery when connected to a temporary external cardiac pacemaker or recording device.

CONTRAINDICATIONS

The device is not intended for contraceptive tubal coagulation (permanent female sterilization).

ENVIRONMENT SPECIFICATIONS

Operational	Storage	Transit
Temperature: 10°C/50°F to 40°C/104°F	Temperature: -29°C/-20°F to 60°C/140°F	Temperature: -29°C/-20°F to 60°C/140°F
Relative Humidity: 15% to 90%, non-condensing	Relative Humidity: 15% to 85%	Relative Humidity: 30% to 85%
Atmospheric Pressure: 80.0 kPa/11.6 psi to 106 kPa/15.4 psi	N/A	N/A

PACKAGING CONTENTS

- One (1) AtriCure Isolator Pen

The Isolator Pen is supplied STERILE and is NON-PYROGENIC in unopened, undamaged package. For single use only. Do not re-sterilize. Do not re-use.

PROVIDED SEPARATELY

- OSCOR PACE 203H and ORLAB
- AtriCure RF Generator

NOMENCLATURE

This instruction refers to features of the Isolator Pen (MAX3) (Figure 1) as follows.

ISOLATOR PEN (MAX3) FEATURES

[1] Distal Tip	[5] Cable
[2] Electrodes	[6] Connector
[3] Malleable Shaft	[7] Alignment Arrow (See Figure 2)
[4] Handle	

POTENTIAL COMPLICATIONS

DEVICE

Possible complications related to the creation of spot or linear lesions in cardiac and soft tissues are:

- Ablation or burns to non-targeted tissues
- Tissue perforation
- Postoperative embolic complications
- Extension of extracorporeal bypass
- Perioperative heart rhythm disturbance (atrial and/or ventricular)
- Pericardial effusion or tamponade
- Damage to adjacent nerve and/or blood vessels
- Valve leaflet damage
- Conduction disturbances (SA/AV node)
- Acute ischemic myocardial event

PROCEDURE

⚠ WARNINGS ⚠

To avoid shock/burn hazards/ineffective cardioversion, always remove the Pen from the patient during defibrillation.

Use of the Pen should be limited to properly trained and qualified medical personnel. Proper surgical procedures and techniques are the responsibility of the medical professional. Understanding the proper use of the Oscor PACE 203 H temporary pacemaker equipment is also the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of any procedure based on their own medical training and experience, and the type of surgical procedure.

Do not use the Pen in the presence of flammable anesthetics, other flammable gases, flammable cleaning agents, near flammable fluids such as skin prepping agents and tinctures, flammable objects, or with oxidizing agents. Use near flammable agents may result in fire or explosion. Observe appropriate fire precautions at all times.

Do not use the Pen for coagulation or ablation of veins or arteries.

Read all instructions carefully for the AtriCure Generator, Isolator Transpolar Pen, Generator Source Switch, and any auxiliary device being used prior to using the devices. Failure to properly follow instructions may lead to electrical or thermal injury and may result in improper functioning of the device.

⚠ CAUTIONS:

- Do not touch the Electrodes of the Pen to metal staples or clips, or to sutures while activating the Generator. This may damage the Pen or tissue or result in an incomplete ablation.
- Do not immerse the Pen in liquids as this may damage the device.
- Do not allow the Connectors of the Pen to get wet. Wet Connectors may affect the device performance.

INSTRUCTIONS FOR USE SET UP

⚠ WARNINGS ⚠

To avoid the risk of patient infection, inspect the product packaging prior to opening to ensure that the sterility barrier is not breached. If the sterility barrier is breached, do not use the Pen.

To avoid damage to the device, do not drop or toss the Pen. If the Pen is dropped, do not use. Replace with a new Pen.

- Inspect the package and product to ensure the expiration date has not passed and no damage occurred to the product during shipping and handling. Using sterile technique, remove the Pen from its packaging.

⚠ WARNINGS ⚠

Use only connecting cables and auxiliary device identified in this Instruction for Use to avoid risk of patient injury, operator injury, or equipment damage.

Do not connect the Generator cable to supply mains (line voltage) operated equipment without verifying isolation of the connected equipment to BS EN 60601-1-1.

Supply mains operated equipment may introduce dangerous leakage currents into the heart.

⚠ CAUTION: The Pen is only compatible with the AtriCure Generator. Use of the Pen with another manufacturer's generator may damage the device and result in patient injury.

⚠ CAUTION: Use care when connecting the Pen. Ensure connector is fully and properly seated in the pen port.

- With the Connector Alignment Arrow symbol in the 12 o'clock position, push the Connector into the Pen receptacle on the front of the Generator (See Figure 2). Verify that the connections between the Pen and the Generator are secure. If the Connections are loose, do not use the Pen. Inspect the Cable and do not use the Pen if the cable is frayed or the insulation is damaged.

PACING AND MAPPING MODE

- Connect one end of the auxiliary device cable provided with the Generator into the compatible auxiliary device. Connect the other end of the auxiliary device cable into one set of PSS ports on the Generator.
- Rotate the Generator knob to indicate the PSS port used.

- Turn on the temporary pacing or sensing equipment and ensure proper connections to validate electrical continuity. For detailed instructions refer to temporary pacemaker manual.
- Settings and procedures for the auxiliary device are determined according to the instructions for use provided with the auxiliary device. Set auxiliary device to atrial asynchronous mode (sensing disabled or increased to maximum value). The Pen will pace when the auxiliary device is in the ON position.

⚠ CAUTION: Ensure the Electrodes are placed on the target tissue. Applying pacing, stimulation, or sensing to unintended structures will give incorrect results.

⚠ CAUTION: Use caution during device insertion and removal.

⚠ CAUTION: Excessive bending of the Malleable stainless-steel Shaft will cause the Shaft to harden and may increase the potential for breakage.

- Identify and expose the sites for pacing and sensing using standard surgical techniques. Under direct visualization, place the Electrodes against the targeted tissue. Assure both Electrodes are in contact with targeted tissue. If necessary, bend the malleable zone of the Shaft to orient the end-effector.
- Ensure the pacing is only on when the end effector is in contact with the targeted tissue. The device will pace when the auxiliary device pacing is on.
- For sensing (mapping), place the Distal Tip onto the targeted tissue to display the electrogram (EGM).
- Upon completion of the surgical procedure, disconnect the Pen from the Generator and discard the Pen after use. See device disposal section below.

ABLATION

- Rotate the Generator knob to indicate the Pen port used.

⚠ WARNINGS ⚠

Ensure constant firm pressure is applied to the tissue without movement.

Application of light pressure, excessive pressure, uneven pressure, or movement of Pen can result in tissue damage or non-transmural lesions.

Use caution to avoid trauma to tissues not within the target area of ablation.

Tissue and/or structures behind the targeted tissue should be protected from potential thermal spread.

Ensure the Electrodes are placed on the target tissue.

Ablation of unintended tissues or structures can cause harm and damage to patient.

⚠ Caution: Do not touch the Electrodes of the Pen while activating the Generator. Touching the Pen Electrodes during Generator activation could result in an electrical shock or burn to the operator.

- Painting ablation technique

- Maintaining visualization, move the Distal Tip gently across the targeted cardiac tissue. If necessary, bend the malleable zone of the Shaft to orient the end-effector.
- While maintaining continuous contact between the tissue and the Electrodes, move the device continuously in an oscillating manner at a rate of approximately 1cm/sec.

Painting Lesion Thickness*
20 seconds
2 cm oscillation
2.0 - 4.0 mm

*Data was obtained from ablations performed on excised bovine myocardium and represent 95% confidence intervals. Results may vary based on live tissue properties.

⚠ WARNINGS ⚠

Ensure stamping lesions are overlapping by 50% to ensure a continuous and complete lesion.

Failure to overlap stamping lesions may result in an undetected, incomplete lesion.

Overlapping by more than 50% may result in tissue perforation or unintended damage.

Total duration of ablation(s) per lesion not to exceed recommended ablation time.

Do not overlap ablations by more than 50% to avoid tissue damage.

- Stamping Ablation technique:

- Apply constant firm pressure to the tissue without movement. Maintain full contact of the Electrode surface with the tissue. If necessary, bend the malleable zone of the Shaft to orient the end-effector. A stamping lesion is approximately 8 mm x 6 mm.
- If creating longer lesions with the Stamp technique, overlap the contiguous ablations by 50% to ensure a continuous and complete lesion.

Stamp Lesion Thickness*	
10 seconds	15 seconds
3.3 - 3.8 mm	3.8 - 4.4 mm

*Data was obtained from ablations performed on excised bovine myocardium and represent 95% confidence intervals. Results may vary based on live tissue properties.

⚠ WARNINGS ⚠

Ensure the Footswitch or button is depressed fully and held down for the entirety of the desired RF energy delivery.

To avoid incomplete lesions of ablation of unintended tissues or structures, use care to depress Footswitch or button only when desired and when the device is placed on target tissue.

Ensure that energy is applied to tissue for the proper amount of time. Energy applied for an extended period of time may cause lateral lesion to spread or tissue damage. Energy applied for a shortened period of time may cause an incomplete lesion.

- Press and hold the Footswitch or button to activate the Generator.
- When the Footswitch or button is pressed, the Generator will emit an audible tone indicating that current is flowing between the Electrodes located at the Distal Tip of the Pen and through the tissue.

⚠ WARNINGS ⚠

Inspect the surgical area to ensure adequate ablation. Failure to inspect may result in an undetected, incomplete lesion.

The Distal Tip of the Pen must be kept clean of debris to avoid loss of power. Before activating the Generator, inspect the area at the Distal Tip of the Pen for foreign matter. Foreign matter captured on the tip will adversely affect the ablation.

⚠ CAUTION: Do not use abrasive cleaners or electrosurgical tip cleaners to clean debris from the Distal Tip. Do not activate the Generator while cleaning the device. Use of abrasive cleaners or electrosurgical tip cleaners can damage the Electrodes and result in device failure. Use saline-soaked gauze for cleaning debris.

- Inspect the surgical area to ensure adequate ablation. Between ablations, wipe the Distal Tip clean with a saline soaked gauze pad. For optimal performance, keep the Pen Electrodes clear of coagulum. The coagulum is much easier to remove within the first several seconds after ablation. In a brief period of time, the coagulum may dry making removal of coagulum more difficult. Check both Electrodes before each ablation to ensure that the gold of each Electrode is visible and coagulum is removed.
- If the Pen is idle between ablations, place the Pen Distal Tip onto saline soaked gauze to prevent any coagulum not cleaned off the Electrodes from drying. Repeat ablation if necessary.

⚠ WARNINGS ⚠

The Pen device is intended for single use only. Do not RESTERILIZE. Re-sterilization may cause loss of function or injury to patient. It is the responsibility of the user to dispose of this device in accordance with local regulations.

⚠ CAUTION: The Pen has an eight-hour useful life that is tracked by the Generator. If you attempt to plug in a device that has reached its time limit expiration, the Pen will no longer function, and the Generator will display a message indicating that the Pen must be replaced.

- Upon completion of the surgical procedure, disconnect the Pen from the Generator and discard the Pen after use. After use, this device should be treated as medical waste and disposed of following hospital procedure. See section below on device disposal.

REMOVAL AND DISPOSAL

After use this device should be treated as medical waste and disposed of following hospital protocol.

RETURN OF USED PRODUCT

If for any reason these products must be returned to ATRICURE, a return goods authorization (RGA) number is required from ATRICURE prior to shipping.

If the products have been in contact with blood or body fluids, they must be thoroughly cleaned and disinfected before packing. They should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and they should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from ATRICURE, Inc.

It is the responsibility of the health care institution to adequately prepare and identify the products for shipment.

DISCLAIMER STATEMENTS

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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