
**Rx ONLY**

# AtriCure®

## ABLATION AND SENSING UNIT (ASU)

### INSTRUCTIONS FOR USE

# ASU2-115 / ASU3-230

### INSTRUCTIONS FOR USE

**MD**

**⚠ Caution:** Federal Law (US) restricts this device to sale by or on the order of a physician.

AtriCure Inc.  
7555 Innovation Way  
Mason, Ohio 45040 USA  
USA  
+1 866 349 2342  
+1 513 755 4100

Commercial Name: AtriCure Synergy Ablation System

Technical Name: RF Ablation Equipment

Contents:

01 Unit IFU

01 Unit of the AtriCure Ablation and Sensing Unit

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Email: [brazilvigilance@ul.com](mailto:brazilvigilance@ul.com)

**THIS PAGE INTENTIONALLY LEFT BLANK**

# TABLE OF CONTENTS

<b>1. Getting Started</b> .....	<b>1</b>
1.1. System Description .....	1
1.2. Indications for Use .....	1
1.3. Unpacking .....	1
1.4. Warnings and Cautions .....	1
1.5. EMC Guidance and Manufacturer's Declaration .....	3
1.6. Responsibility of the Manufacturer .....	6
<b>2. The AtriCure Ablation and Sensing Unit (ASU)</b> .....	<b>6</b>
2.1. Device Description .....	6
2.2. ASU Front Panel – Illustration and Nomenclature .....	6
<b>Front Panel Displays</b> .....	<b>6</b>
<b>Front Panel Indicators</b> .....	<b>7</b>
<b>Front Panel Receptacle</b> .....	<b>7</b>
2.3. ASU Rear Panel – Illustration and Nomenclature.....	8
<b>Rear Panel Functions</b> .....	<b>8</b>
<b>3. Installing the ASU</b> .....	<b>9</b>
3.1. Transporting the ASU.....	9
3.2. Adjusting the Viewing Angle .....	9
3.3. Preparing the ASU For Use .....	9
3.4. Power Cord .....	9
3.5. Connecting and Disconnecting the Handpiece .....	9
3.6. Installing the Footswitch .....	10
<b>4. Instructions For Use</b> .....	<b>10</b>
4.1. Powering Up the ASU .....	10
4.2. Operating Modes .....	11
4.3. Audio Tones .....	11
4.4. Delivering RF Energy .....	12
<b>5. Troubleshooting</b> .....	<b>14</b>
5.1. No RF Power Output .....	14
5.2. Error Codes .....	14
5.3. Electromagnetic or Other Interference .....	15
<b>6. Symbols Used</b> .....	<b>16</b>
<b>7. TECHNICAL SPECIFICATIONS AND SAFETY INSPECTION</b> .....	<b>16</b>
7.1. RF Output .....	16
7.2. Mechanical Specifications .....	17
7.3. TECHNICAL SPECIFICATIONS .....	17
7.4. Environmental Specifications .....	17
7.5. Electrical Specifications .....	17
7.6. Fuse Specifications .....	17
7.7. Footswitch Specifications .....	17
7.8. Power and Voltage Output Restrictions .....	17
7.9. Equipment Type / Classification .....	18
<b>8. Preventive Maintenance and Cleaning of ASU</b> .....	<b>18</b>
8.1. Preventive Maintenance .....	18

8.2. Cleaning and Disinfecting.....	19
Guidelines .....	19
<b>9. Disposal .....</b>	<b>19</b>
<b>10. Accessories .....</b>	<b>20</b>
10.1. ASB3, Switch Matrix Accessory .....	20
<b>11. Components and Cables .....</b>	<b>21</b>
Limitation on Liability .....	22
<b>DISCLAIMER .....</b>	<b>22</b>

## 1. GETTING STARTED

This manual, and the equipment, it describes are for use only by qualified medical professionals trained in the particular technique and surgical procedure to be performed.

Please read all information carefully. Failure to properly follow the instructions may lead to serious surgical consequences.

**Important:** This manual is designed to provide instructions for use of the AtriCure Ablation and Sensing Unit (ASU) with the AtriCure Bipolar Handpiece (Synergy™ clamp, Isolator™ Transpolar™ pen) and AtriCure Accessory Devices (ASB3). It is not a reference to surgical technique.

The AtriCure® ASU produces and delivers RF energy, in a bipolar mode, at a frequency of approximately 460 kHz, with a maximum output power ranging from 22.8 Watts up to 28.5 Watts for the Synergy Clamps, 12.0 Watts up to 30.0 Watts for the Isolator™ Transpolar™ pen devices depending on the mode of operation. The AtriCure® ASU is capable of producing a maximum output power of 32.5 Watts under a 100 Ohm load, although no current AtriCure® Bipolar Handpiece uses power above 30 Watts. The operating mode is a function of the handpieces or pen and is set by the ASU. The AtriCure ASU is designed to operate only with an AtriCure Bipolar Handpiece, or AtriCure Isolator Pen. The Footswitch is the input device used to activate RF energy delivery. Please refer to the handpiece and pen Instructions for Use for complete description of the indications and use of these devices.

For the user's convenience, the AtriCure Ablation and Sensing Unit will be referred to in this User's Manual as the "ASU". The AtriCure Bipolar Handpiece will be referred in this User's Manual as the "Handpiece".

This User's Manual provides a description of the ASU, its controls, displays, indicators, tones and a sequence for its operation with the Handpiece. This User's Manual also supplies other information of importance to the user. This manual is intended as a User's Manual only. Do not operate the ASU before thoroughly reading this manual.

### 1.1. System Description

As shown in Figure 1, the system is comprised of the following:

- AtriCure Bipolar Handpiece with integral cable (not shown)
- AtriCure Ablation and Sensing Unit (ASU)
- Footswitch
- Power cord.

Accessory devices are described in paragraph 10.


Figure 1 – ASU, Footswitch, and Power Cord

### 1.2. Indications for Use

The AtriCure Ablation and Sensing Unit (ASU) is designed to provide energy to various AtriCure RF ablation handpieces (pens and clamps) that are intended for surgical ablation of cardiac tissue.

### 1.3. Unpacking

Lift the ASU, Footswitch, and Power Cord from the box and remove the protective wrapping. It is recommended that the original shipping box and protective wrapping be saved for future storing and/or transporting of the device.

### 1.4. Warnings and Cautions

The safe and effective use of RF energy is highly dependent upon factors under the control of the operator. There is no substitute for a properly trained operating room staff. It is important that the operating instructions supplied with the ASU be read, understood and followed before use.

#### WARNINGS

- Do not operate the ASU/ASB before thoroughly reading this manual.
- Do not use the ASU/ASB in the presence of flammable anesthetics, other flammable gases, flammable cleaning agents, near flammable fluids such as skin prepping agents and tinctures, flammable objects, or with oxidizing agents. Use near flammable agents may result in fire or explosion. Observe appropriate fire precautions at all times.

- Fire Hazard: Electrosurgical accessories that are activated or hot from use can cause a fire. Do not place them near or in contact with flammable materials (such as gauze or surgical drapes). Avoid igniting endogenous gases.
- Fire Hazard: To avoid igniting cleaning agents, use only non-flammable agents to clean and disinfect the ASU/ASB. If flammable agents are inadvertently used on the Generator, allow these substances to evaporate completely before operating.
- Electrosurgery should be used with caution in the presence of internal or external pacemakers. Interference produced with the use of electrosurgical devices can cause devices such as a pacemaker to enter an asynchronous mode or can block the pacemaker entirely. Consult the pacemaker manufacturer or hospital Cardiology department for further information when the use of electrosurgical appliances is planned in patients with cardiac pacemakers.
- The use of components, transducers, and cables other than those specified in accordance with the instructions or supplied by AtriCure, may result in increased emissions or decreased immunity of the equipment.
- The ASU/ASB should not be used adjacent or stacked with other equipment, except for intended stacking with AtriCure's equipment in accordance with the instructions. The ASU/ASB normal use configuration should be observed to verify normal operation.
- The voltage selector is factory set and should not be changed by the user. The voltage selector and the power entry module must be set to the same voltage setting to prevent ASU malfunction and potential instrument damage.
- The Power cord(s) of the ASU/ASB must be connected to a properly grounded receptacle(s). Extension cords and/or adapter plugs must not be used.
- Electric Shock Hazard: Do not connect wet components to the Generator.
- Electric Shock Hazard: Ensure that the Handpiece is correctly connected to the Generator and that no wires are exposed from the cable, connector, or Handpiece.
- Do not connect the ASU/ASB device cable to supply mains (line voltage) operated equipment without evidence that the safety certification of the Accessory has been performed in accordance to the appropriate EN60601-1 and/or EN60601-1-1 harmonized national standard. Supply mains operated equipment may introduce dangerous leakage currents into the heart. When transporting or handling the ASU/ASB, use caution and care to avoid product damage. Inspect the ASU/ASB and Interface cables for any signs of physical damage before use. If damage is found, the integrity of the product cannot be guaranteed, and the product should not be used.
- Make sure there are no obstructions underneath or near the rear of the ASU/ASB to provide sufficient air flow for cooling.
- Use care when connecting Footswitch, Power cords, or Handpieces. Failure to connect properly may result in product damage.
- The safe and effective use of RF energy is highly dependent upon factors under the control of the operator. The manual, and the equipment it describes, are for use only by qualified medical professionals trained in the particular technique and surgical procedure to be performed. It is important that the operating instructions supplied with the ASU/ASB be read, understood and followed before use to avoid product damage or patient harm.
- If the ASU powers down, use SECTION 5 TROUBLESHOOTING to attempt to correct the problem. If problem persists, contact the AtriCure representative.
- The activation tone and indicator are important safety features. Do not obstruct the activation indicator or disable the audible tone. Ensure that the activation tone is audible to personnel in the Operating Room (OR) prior to use. The activation tone alerts personnel when the Handpiece is active and can help avoid tissue damage.
- Follow the cleaning instructions detailed in SECTION 8: PREVENTIVE MAINTENANCE AND CLEANING. Failure to follow these guidelines may result in product damage.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ASU/ASB including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

## Cautions

- Use only with the AtriCure Handpieces intended for use with the ASU/ASB.
- Do not activate the Generator until the Handpiece is properly positioned in the patient.
- Electric shock hazard: Do not remove the cover of the ASU/ASB as there is a potential for electric shock. Refer to authorized personnel for service.
- Trip Hazard: Standard care should be used to reduce the risk of tripping on the Footswitch cable.
- Use only the Footswitch provided with the ASU/ASB.
- Do not wrap instrument cables around metal objects. Wrapping cables around metal objects may induce hazardous currents.
- To avoid shock, do not allow patients to come into contact with earth metal parts of the Device. The use of antistatic sheeting is recommended.
- When the ASU/ASB is activated, the conducted and radiated electrical fields may interfere with other electrical medical equipment. Refer to the EMC Guidance and Manufacturer's Declaration section for more information regarding potential electromagnetic or other interference and advice regarding avoidance of such interference.
- Studies have shown that smoke generated during electrosurgical procedures can be potentially harmful to surgical personnel. These studies recommend using surgical masks and adequately ventilating the smoke by using a smoke evacuator or other means.
- When the Generator and Handpiece are used on a patient simultaneously with physiological monitoring equipment, ensure that the monitoring electrodes are placed as far as possible from the surgical electrodes. Be sure to position the Handpiece cables so that they do not come in contact with the patient or the other leads.
- Needle monitoring electrodes are not recommended for use when operating the Generator and the Handpiece.
- Monitoring systems that incorporate high-frequency current-limiting devices are recommended for use with the Generator and the Handpiece.
- Failure of the Device and the Handpiece could result in unintended Power output increases.
- Handle the ASU and ASB packaging with care.
- The ASU/ASB, AtriCure Isolator Transpolar Pen, Isolator Synergy Clamps have been tested as a system. Another manufacturer's accessories are not compatible and will not connect to the ASU/ASB.
- The ASU/ASB generates, uses, and can radiate RF energy. Interference produced by the operation of the ASU/ASB may adversely influence the operation of other electronic medical equipment such as monitors and imaging systems.

## 1.5. EMC Guidance and Manufacturer's Declaration

### 1.5.1. Electromagnetic Requirements

The AtriCure Ablation and Sensing Unit (ASU) has been tested and found to comply with the limits for medical devices in EN 60601-1-2:2015 + A1:2021. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This system generates, uses, and can radiate radiofrequency energy and, if not installed and used in accordance with the instruction given below, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. The ASU/ASB do not require preventive maintenance with regard to electromagnetic disturbance during the expected service life. If the ASU/ASB does cause harmful interference to other devices, which can be determined by turning the ASU/ASB OFF and ON, the operator is encouraged to try to correct the interference by:

- Relocate or reorient the equipment
- Increase the separation distance between the equipment and the other devices
- Connect the equipment into other outlets different from those to which the other devices are connected
- Consult AtriCure, Inc. representatives for help.

### 1.5.2. Electromagnetic Emissions

Table A: IEC EMC Specifications (Emissions)		
<b>Guidance and Manufacturer's Declaration – Electromagnetic Emissions</b>		
The AtriCure Ablation and Sensing Unit (ASU) is intended for use in the electromagnetic environment specified below. The customer or the user of the ASU System should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
Radiated RF emissions CISPR 11	Group 1 Class A	The AtriCure Ablation and Sensing Unit (ASU) must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
Conducted RF emissions CISPR 11	Group 1 Class A	The AtriCure Ablation and Sensing Unit (ASU) is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics emissions IEC 61000-3-2	Class A	
Voltage changes / fluctuations / flicker IEC 61000-3-3	Complies	

**NOTE:** The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required), this equipment might not offer adequate protection to radiofrequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

### 1.5.3. Electromagnetic Immunity

Table B: IEC EMC Specifications (Immunity)			
<b>Guidance and Manufacturer's Declaration – Electromagnetic Immunity</b>			
The AtriCure Ablation and Sensing Unit (ASU) is intended for use in the electromagnetic environment specified below. The customer or the user of the ASU System should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV @ 100 kHz repetition frequency for power supply lines ± 2 kV @ 100 kHz repetition frequency for input/output lines	± 2 kV @ 100 kHz repetition frequency for power supply lines ± 2 kV @ 100 kHz repetition frequency for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Power inputs ± 0,5 kV, ± 1 kV Line-to-Line ± 0,5 kV, ± 1 kV, ± 2 kV Line-to-Ground Signal input/outputs: ± 2 kV Line-to-Ground	Power inputs ± 0,5 kV, ± 1 kV Line-to-Line ± 0,5 kV, ± 1 kV, ± 2 kV Line-to-Ground Signal input/outputs: ± 2 kV Line-to-Ground	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage Dips: 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° Voltage interruptions: 0 % UT; 250/300 cycle	Voltage Dips: 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0 % UT; 1 cycle Single phase: at 0° 70 % UT; 25/30 cycles Single phase: at 0° Voltage interruptions: 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ASU System requires continued operation during power mains interruptions, it is recommended that the ASU System be powered from an uninterruptible power supply or battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	0,15 MHz – 80 MHz 3V, 80 % AM at 1 kHz  ISM bands between 0,15 MHz and 80 MHz 6V, 80 % AM at 1 kHz	0,15 MHz – 80 MHz 3V, 80 % AM at 1 kHz  ISM bands between 0,15 MHz and 80 MHz 6V, 80 % AM at 1 kHz	Mains power quality should be that of a typical commercial or hospital environment.

**NOTE** UT is the a.c. mains voltage prior to application of the test level.

#### 1.5.4. EMC Guidance and Manufacturer's Declaration

**Table C: IEC EMC Specifications (Immunity from Radiated RF EM Fields)**

Immunity test	Band (MHz)	Wireless Service	Immunity Test Level (V/m)	Compliance Test Level (V/m)
Immunity from Radiated RF EM Fields including proximity fields from RF wireless communications equipment  IEC 61000-4-3	150 kHz to 80 MHz	General	< 3	< 3
	80 MHz – 2,7 GHz	General	3	3
	380 – 390	TETRA 400	27	27
	430 – 470	GMRS 460, FRS 460	28	28
	704 – 787	LTE Band 13, 1	9	9
	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	28	28
	1,700 – 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	28	28
	2,400 – 2,570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	28	28
5,100 – 5,800	WLAN 802.11 a/n	9	9	

**Table C: IEC EMC Specifications (Immunity from Radiated RF EM Fields)**

Immunity test	Band (MHz)	Wireless Service	Immunity Test Level (V/m)	Compliance Test Level (V/m)
<p>Portable and mobile RF communications equipment should be used no closer to any part of the ASU System including cables, than the recommended separation distance calculated from the equation:</p> <p style="text-align: center;">Where:</p> $d = 6/E \times \sqrt{P}$ <p style="text-align: center;"><i>d</i> is the separation in meters</p> <p style="text-align: center;"><i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the Service</p> <p style="text-align: center;"><i>E</i> is the Compliance Test Level indicated above.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  <p>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ASU System or any of its components are used exceeds the applicable RF compliance level above, the ASU System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating components or the entire ASU System.</p> <p style="text-align: center;">Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>				

**Table D: Immunity to Proximity Magnetic Fields**

Test Frequency	Modulation	Immunity Test Level (A/m)
30 kHz <sup>a</sup>	CW	8
134,2 kHz	Pulse modulation <sup>b</sup> 2,1 kHz	65 <sup>c</sup>
13,56 MHz	Pulse modulation <sup>b</sup> 50 kHz	7,5 <sup>c</sup>

a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT.  
 b) The carrier shall be modulated using a 50% duty cycle square wave signal.  
 c) r.m.s, before modulation is applied.

**1.5.5. Recommended Separation Distance**

**Recommended separation distances between portable and mobile RF communications equipment and the AtriCure Ablation and Sensing Unit**

The AtriCure Ablation and Sensing Unit (ASU) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ASU can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ASU as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## 1.6. Responsibility of the Manufacturer

AtriCure is responsible for safety, reliability, and performance of the equipment only if:

- Installation procedures in this manual are followed.
- Persons authorized by AtriCure carry out modifications or repairs.
- The electrical installation of the relevant room complies with local codes and regulatory requirements such as IEC and BSI.
- The equipment is used in accordance with the AtriCure User's Manual.

## 2. THE ATRICURE ABLATION AND SENSING UNIT (ASU)

This section provides a detailed description of the ASU including its function and operating features.

### 2.1. Device Description

The AtriCure® ASU produces and delivers RF energy, in a bipolar mode, at a frequency of approximately 460 kHz, with a maximum output power ranging from 12 Watts up to 30 Watts depending on the operating mode. The AtriCure® ASU is capable of producing a maximum output power of 32.5 Watts under a 100 Ohm load although no current AtriCure® Bipolar Handpiece uses power above 30 Watts. The operating mode is a function of the handpiece and is set by the ASU. The AtriCure ASU is designed to operate with the AtriCure Handpiece. The ASU and Handpiece are designed for use without a neutral electrode. The Footswitch is the input device used to activate RF energy delivery.

### 2.2. ASU Front Panel – Illustration and Nomenclature

An illustration of the ASU front panel is shown in Figure 2, below.



Figure 2 – ASU Front Panel

- |  |                             |
|--|-----------------------------|
| 1. Tissue Conductance/ Power Graph Display | 6. Transmurality Indicator  |
| 2. Ready Indicator                         | 7. Fault Indicator          |
| 3. RF ON Indicator                         | 8. Power Indicator          |
| 4. Temperature Display                     | 9. Handle                   |
| 5. Handpiece Receptacle                    | 10. Handle Adjustment Knobs |

### Front Panel Displays

There are two displays on the front panel of the ASU: the Tissue Conductance / Power Graph Display and the Temperature Display. These two displays are described below.


Display	Description
	<p><b>Tissue Conductance Graph Display – Isolator® Synergy™ clamp (Default):</b></p> <p>During the ablation cycle the ASU displays a graph of tissue conductance (Current/Voltage) versus Time. The y-axis is Tissue Conductance and the x-axis is Time.</p> <p>When the Footswitch is disconnected or reconnected, the display of the tissue conductance graph is not affected. Refer to Section 4.4.3.</p>

Display	Description
	<p><b>Power Graph Display – Isolator™ Transpolar™ pen:</b></p> <p>During the ablation cycle the ASU displays a graph of power (Current™ Voltage) versus Time. The y-axis is Power and the x-axis is Time.</p> <p>When the Footswitch is disconnected or reconnected, the display of the power graph is not affected. Refer to Section 4.4.3.</p>
	<p><b>Temperature Display</b> – This 3-digit LED display shows the temperature at the thermocouple, located near the outer edge of the upper jaw 1.3 mm from the electrode. Temperature is measured and displayed in real time, whenever the Isolator® Synergy™ clamp is connected. Functionality can be quickly verified when the handpiece is plugged in by confirming that the temperature reading is of room temperature.</p> <p>When the Isolator® Synergy™ clamp or the footswitch is disconnected the temperature display becomes blank. Refer to Section 4.3.</p> <p>If a Handpieces does not have a thermocouple, the temperature display will only show “- - -”.</p>

### Front Panel Indicators

Indicator	Description
	POWER Indicator – A Green LED indicates that the AC power is present and the ASU has been switched on.
	FAULT Indicator – This Red lamp indicates that a fault has occurred and requires that the power be cycled.
<p>READY</p> 	READY Indicator – This Green lamp indicates that the Footswitch and Handpiece are connected and the ASU is ready for use
<p>RF ON</p> 	RF ON Indicator – A Blue LED indicates that RF power is being output to the Handpiece.
<p>TRANSMURALITY</p> 	TRANSMURALITY Indicator – A Blue flashing LED indicates that the Transmurality Algorithm has been satisfied indicating that the user may terminate the ablation cycle.

### Front Panel Receptacle

Receptacle	Description
	<p><b>HANDPIECE or ASU Accessory Receptacle</b></p> <p>This 12-pin receptacle accepts the AtriCure Handpiece or connection cable to an accessory device. This connection is patient-isolated.</p>

## 2.3. ASU Rear Panel – Illustration and Nomenclature

An illustration of the ASU rear panel is shown in Figure 3, below.



FIGURE 3 – ASU REAR PANEL

- |                           |                                  |
|---------------------------|----------------------------------|
| 1. Data Port              | 5. Equipotential Ground Stud     |
| 2. Speaker Volume Control | 6. Input Voltage Selector Switch |
| 3. Power Entry Module     | 7. Footswitch Receptacle         |
| 4. Fuse Box               | 8. Service Access                |

### Rear Panel Functions

Graphic	Description
	Equipotential Ground Stud – Provides a means of securely linking the earth grounds of the AtriCure ASU to other grounded equipment.
	Data Port – For manufacturing and test purposes.
	Power Entry Module – This module contains both the ON/OFF switch and the fuses. The voltage is selected by the orientation of the fuse drawer as marked.  Fuse Box – The Fuse Box contains fuses selected for the input voltage. See Technical Specifications in Section 7 of this manual.
	Input Voltage Selector Switch – The input voltage selector switch is pre-set at the factory to either 110V or 220V and should not be adjusted by the operator. This setting should only be adjusted by the manufacturer or by an authorized service representative.

Graphic	Description
<p><b>VOLUME</b></p> 	<p>Speaker Volume Control – The audible volume level is adjustable via a rotary dial.</p> <p>The ASU includes a speaker for producing audible feedback to the user.</p>
<p><b>FOOTSWITCH</b></p> 	<p>Footswitch Receptacle – This receptacle accepts the Footswitch connector. The single momentary actuation pedal provides for the activation of RF power output.</p>
<p><b>SW1</b></p> 	<p>Service Access – For manufacturing and test purposes.</p>

### 3. INSTALLING THE ASU

Use Gloves when setting up and operating the ASU/ASB system and associated components.

Inspect the ASU for any signs of physical damage to the front panel, chassis or cover.

**NOTE:** If any physical damage is found, DO NOT USE THE UNIT. CONTACT AtrıCure for a replacement.

All returns must be approved by AtrıCure.

#### 3.1. Transporting the ASU

The handle may be used to carry the ASU. To change the positioning of the handle, depress both handle adjustment knobs simultaneously and move the handle to the desired location. **Do not** change the handle position when a Handpiece or Accessory Device is connected to the Handpiece receptacle.

#### 3.2. Adjusting the Viewing Angle

To change the viewing angle of the ASU Conductance Graph Display, adjust the handle position using the directions in Section 3.1., above.

#### 3.3. Preparing the ASU For Use

The ASU may be placed on a mounting cart or on any sturdy table or platform. It is recommended that carts have conductive wheels. Refer to hospital procedures or local codes for detailed information.

Provide at least four to six inches of space around the sides and top of the ASU for convection cooling. Under continuous use for extended periods of time, it is normal for the top and rear panel to be warm.

Ensure there are no obstructions underneath or near the back of the ASU/ASB to provide sufficient airflow for cooling. Connect cables and handpieces into receptacles with care to avoid damage.

#### 3.4. Power Cord

The ASU is shipped with an approved hospital grade power cord.

 **Caution:** Inspect cord insulation, connector ends, and receptacles for damage before use.

Plug the ASU into a grounded receptacle.

Do not use extension cords or three-prong to two-prong adapters. The Power Cord assembly should be periodically inspected for damages per SECTION 8 PREVENTIVE MAINTENANCE AND CLEANING OF ASU.

#### 3.5. Connecting and Disconnecting the Handpiece

Connect the Handpiece directly to the ASU. Insert the Handpiece cable connector into the receptacle on the front panel of the ASU, ensuring that the arrow symbol on the connector is facing upward and oriented to the arrow symbol on the ASU receptacle.

**NOTE:** Typically, you will connect the Handpiece to the ASU when the ASU has been powered up and is in STANDBY operating mode (see Section 4.2 regarding the STANDBY mode). However, the Handpiece may be connected when powered up, or prior to powering up the ASU.

**NOTE:** Once you have connected the Handpiece, it cannot be disconnected from the ASU by pulling on the cable. To disconnect the Handpiece, pull back on the cable connector body and remove it from the ASU receptacle.

**NOTE:** Refer to the Handpiece instruction sheet for more detailed information about connecting the Handpiece to the ASU in a sterile environment.

### 3.6. Installing the Footswitch

#### 3.6.1. Inspect the Footswitch

Inspect the Footswitch for any signs of physical damage to the cable and connector. If physical damage is found or the Footswitch does not perform within specification, notify AtriCure. All returns must have approval from AtriCure.

#### 3.6.2. Connecting and Disconnecting the Footswitch

With the connector alignment arrow in the 12 o'clock position, push the Footswitch Connector into the Footswitch Receptacle on the rear panel of the ASU, shown in the Figure 4.

**NOTE:** Typically, you will connect the Footswitch to the ASU when the ASU has been powered up and is in STANDBY operating mode (see Section 4 regarding the STANDBY mode). However, the Footswitch may be connected when powered up, or prior to powering up the ASU.


Figure 4 – Connecting the Footswitch to the ASU

#### 3.6.3. Preparing the Footswitch for Use

The Footswitch should be placed on a flat floor. It is recommended that the area near the Footswitch be kept dry to reduce the risk of slippage.

Appropriate precautions should be taken to ensure that the cable connecting the Footswitch to the ASU does not create a hazard in the operating room.

## 4. INSTRUCTIONS FOR USE

**⚠ Caution:** Ensure the power cord, Footswitch, and Handpieces are installed into the correct receptacle.

### 4.1. Powering Up the ASU

1. Ensure that the ASU has been plugged into a grounded receptacle.

**NOTE:** Do not use extension cords or three-prong to two-prong adapters. The power cord assembly should be periodically checked for damaged insulation or connectors.

2. Turn the power on using the ON/OFF switch located on the power entry module on the rear panel. When power is turned on; the system performs the System Self-Tests. See Figure 5. If all Self-Tests pass, the system transitions to the STANDBY mode. If any Self Test fails, the system transitions to the FAULT mode. The Self-Test generates two quick beeps at startup. The operator must verify that the beeps are generated.

**NOTE:** Refer to Section 4.2., below, for a full description of the STANDBY and FAULT modes, as well as all the other operating modes.

Figure 5 – Display Indicating SELF-TEST

## 4.2. Operating Modes

The ASU operates in one of five modes: STANDBY, READY, RF ON, ERROR and FAULT modes. These modes are shown on the lower left corner of the Conductance Display Graph. See Figure 6, below.

- **STANDBY Mode** – This mode is entered automatically after the ASU is successfully turned on or from READY Mode upon detection of a Handpiece or Footswitch disconnection. The LCD display message indicates the system is in the STANDBY Mode.
- **READY Mode** – This Mode is entered upon connecting both Handpiece and Footswitch while in the STANDBY Mode or from the ON Mode if the Footswitch has been depressed and released. The LCD display message indicates the system is in the READY Mode.
- **RF ON Mode** – This Mode is entered when the Footswitch is depressed while in the READY Mode. The system transitions from the RF ON Mode to the READY Mode upon 40-second time expiration or if the Footswitch is released.
- **ERROR Mode** – This Mode is entered upon detection of any recoverable error conditions during any Mode excluding the FAULT Mode. The system displays the corresponding error message, and upon Footswitch release, transitions to the READY Mode.
- **FAULT Mode** – This Mode is entered upon detection of any unrecoverable error condition during any Mode. The system is inoperable in this Mode until the power is cycled off, then on.


Figure 6 – Conductance and Power Display Graph Indicating STANDBY Mode

**⚠ Caution:** If a maintenance message is displayed on the ASU display screen, contact the AtriCure representative for assistance.

## 4.3. Audio Tones

The ASU uses 7 possible audio tones during its operation: Start Tone, Error Tone, Fault Tone, RF ON Tone, Transmurality Tone, High Temperature RF ON Tone, and the High Temperature Transmurality Tone. You may control the volume of these tones using the Speaker Volume Control on the rear panel of the ASU (See Figure 3). Each of these 7 audio tones is described below.

Tone Name	Tone Description	Meaning for Operator:
Start Tone	Two quick beeps	This tone is generated when the power switch is placed in the “ON” position.

Tone Name	Tone Description	Meaning for Operator:
Error Tone	Constant low-pitched tone	This tone occurs while an error is present.
Fault Tone	Rapid succession of low-pitched beeps for 2 seconds duration	This tone occurs upon entering a fault mode.
RF ON Tone	Constant medium-pitched tone	This tone is generated when RF energy is being delivered to the Isolator® Synergy™ clamp. This tone has a higher pitch than the Error tone.
	Varying medium-pitched tone	A discrete, decrementing tone in 10 second intervals is generated when RF energy is being delivered to the Isolator™ Transpolar™ pen. This tone has a higher pitch than the Error tone.
Transmurality Tone	Intermittent medium-pitched tone	This tone is generated in the RF ON mode when Transmurality is achieved. The Transmurality tone will continue, and RF energy will continue to be applied, until the Footswitch is released or until 40 seconds has elapsed. <b>This function is not applicable to the Isolator™ Transpolar™ pen.</b>
High Temperature RF ON	Constant high-pitched tone	This tone is generated while in the RF ON State when a temperature of 70°C (158°F) or higher is measured by the temperature sensor. The RF energy output is disabled when temperature of 75°C (167°F) or higher is measured by the thermocouple. <b>This function is not applicable to the Isolator™ Transpolar™ pen and the dual electrode Isolator® Synergy™ clamps.</b>
High Temperature Transmurality	Intermittent high-pitched tone	This tone is generated when Transmurality is achieved and when a temperature between 70°C (158°F) to less than 75°C (167°F) is measured by the thermocouple. <b>This function is not applicable to the Isolator™ Transpolar™ pen and the dual electrode Isolator® Synergy™ clamps.</b>

## 4.4. Delivering RF Energy

### 4.4.1. Connect the Handpiece and Footswitch

Connect the Handpiece and Footswitch as described in Sections 3.5. and 3.6., and note the display to ensure connections are made. The display screen and Ready Indicator of the ASU should indicate that the RF generator is in the READY mode. See Figure 7.


Figure 7 – Tissue Conductance Display Graph Indicating READY Mode for Isolator Handpieces(above) and Power Display Graph for Isolator™ Transpolar™ pen device (below).

**NOTE:** When the READY mode is entered from the RF ON mode, the previous plot is shown.

### 4.4.2. Position the Handpiece

To position the Handpiece, follow the Instructions for Use provided with the Handpiece.

**WARNING**

Care should be taken to ensure that the footswitch is not unintentionally activated to avoid ablation of unintended tissue or structures.

**⚠ Caution:** Position the cables to surgical electrodes to prevent contact with the patient or other leads.

### 4.4.3. Deliver RF Energy

**⚠ Caution:** Ensure activation footswitch connected to the ASU/ASB is held down only when delivery of RF energy is desired. Release the footswitch after ablation is completed.

Press the Footswitch to initiate RF energy output. RF energy output is terminated by releasing the Footswitch or at the end of 40 continuous seconds of energy delivery. The display screen of the ASU will indicate that the generator is in the RF ON mode. See Figures 8 and 9.

During the Isolator<sup>®</sup> Synergy<sup>™</sup> clamp operation, a real-time graph of measured tissue conductance is displayed on the LCD graphics screen with a +/- 20% tolerance, and measured temperature is shown on a numeric display. Using measurements of conductance, the ASU will determine when a transmural condition has been achieved.

When this condition has been achieved, the Blue Transmurality indicator will flash and the audible tone emitted from the ASU will change from constant to intermittent, thus signaling to you that transmural condition has been achieved. If you do not release the Footswitch within 40 seconds, the system will automatically time-out and stop the ablation.



Figure 8 – Conductance Display Graph Indicating RF ON Mode

With the Isolator<sup>™</sup> Transpolar<sup>™</sup> pen, a real-time graph of measured power delivered to the tissue is displayed on the LCD graphics screen with a +/- 20% tolerance. The ASU will not indicate when a transmural condition has been achieved in this mode. Furthermore, if you do not release the Footswitch within 40 seconds, the system will automatically time-out and stop the ablation.


Figure 9 – Power Display Graph Indicating RF ON Mode

Both the conductance and the power graphs are on a 20-second scale. In some cases, the transmural condition will not be achieved within the 20 seconds shown on the Tissue Conductance Display Graph (not valid for Isolator<sup>™</sup> Transpolar<sup>™</sup> pen device). In such cases, the Graph will wrap to a second screen, which will display a continuation of the conductance for a maximum of 20 additional seconds. Figure 10, below, shows an example of this wrapping feature for an ablation requiring more than 20 seconds.

Similarly, for the Isolator<sup>™</sup> Transpolar<sup>™</sup> pen the power graph will wrap to a second screen for ablations lasting longer than 20 seconds for a maximum of 20 additional seconds.


Figure 10 – Display Graph Wraps For An Ablation Lasting More Than 20 Seconds

## 5. TROUBLESHOOTING

Use the following sections to help troubleshoot possible problems with the ASU.

**⚠ Caution:** Do not open the back panel of the ASU/ASB. This may cause injury and damage to the units. It will void the warranty. When problems cannot be resolved by the directions in this Troubleshooting section, contact AtriCure, Inc. for additional service and repair information.

### 5.1. No RF Power Output

If there is no RF power output, attempt to correct this problem using the checklist below.

Possible Cause	Solution
ASU not turned on	Turn power on
ASU not plugged in	Confirm electrical connections and then turn power on
No Handpiece connected	Connect Handpiece
No Footswitch connected	Connect Footswitch
ASU in FAULT mode	Turn Power off and then on
ASU in STANDBY mode	Ensure that Handpiece and Footswitch are properly connected
Broken Handpiece cable	Replace Handpiece
Fault in Footswitch	Replace Footswitch
Fault in Handpiece	Replace Handpiece
Internal ASU failure	Contact AtriCure Customer Service

If the lack of ASU RF power output persists, contact the AtriCure service representative.

### 5.2. Error Codes

If a fault condition should occur, the numeric displays on the front panel will display an error code. If an Error Code of E07 through E09, P01 through P11, or F01 through F14 appears, try turning power off, then on. If the problem persists, contact AtriCure Customer Service.

Use the table below to attempt to resolve the following recoverable application errors.

LCD DISPLAY MESSAGE	DESCRIPTION	SOLUTION
Replace Handpiece H01	Invalid Handpiece Version	Replace Handpiece
Replace Handpiece H02	Time Expired Error: The Handpiece expiration date has been exceeded	Replace Handpiece
Replace Handpiece H03	Handpiece Electrical Problem	Replace Handpiece
Replace Handpiece H04	Invalid Handpiece Version	Replace Handpiece
Check Electrodes E01	Low Impedance Error: Handpiece electrodes are shorted	Check Electrodes or reposition jaws
Close Jaws E02	High Impedance Error: Handpiece jaws are open	Close Handpiece Jaws

LCD DISPLAY MESSAGE	DESCRIPTION	SOLUTION
Check Electrodes E03	Low Impedance Error: Handpiece electrodes are shorted	Check Electrodes or reposition jaws
Check Electrodes E04	Low Impedance Error: Handpiece electrodes are shorted	Check Electrodes or reposition jaws
Replace Handpiece E05	Open or defective thermocouple	Replace Handpiece
Check Footswitch E06	Switch Stuck Test Error: Footswitch closed while connecting	Replace Footswitch
Check Electrodes E10	Handpiece electrodes are shorted	Check electrodes or reposition jaws
Return ASU for Maintenance	ASU Clock Battery had failed	ASU will continue to operate, but the message and tone will repeat. Return ASU for battery replacement.

POWER ON SELF TEST (POST) Errors detected at start up		
Code P01	Power Generation / Measurement Error	<p>Power the unit OFF, then Power it ON again.</p> <p>Allow the Device to run through normal start-up self-diagnostics.</p> <p>If the device returns to the Error state and the problem persists, contact AtriCure Customer Service (refer to SECTION 8 CUSTOMER SERVICE/ EQUIPMENT SERVICING/WARRANTY)</p> <p><b>**For P10:</b> Deactivate the footswitch and power cycle the Device.</p>
Code P02	Impedance Generation / Measurement Error	
Code P03	Voltage Generation / Measurement Error	
Code P04	Current Generation / Measurement Error	
Code P05	Watchdog Test Error	
Code P06	ROM Test Error	
Code P07	RAM Test Error	
Code P08	Configuration Register Error	
Code P09	MCU COP Timer Error	
Code P10**	Switch Stuck Test Error (Check Footswitch as it was activated during POST)	
Code P12	Reference Voltage Error	

### 5.3. Electromagnetic or Other Interference

The ASU has been tested and found to comply with the limits for medical devices in EN 60601-1-2:2015 + A1:2021. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

The ASU generates and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If the ASU does cause harmful interference to other devices, which can be determined by turning the ASU off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the ASU and the other devices.
- Connect the ASU into an outlet on a circuit different from that to which the other device(s) are connected.
- Contact the AtriCure service representative for help.

Use the following sections to troubleshoot specific types of interference, including monitor (display) interference, neuromuscular stimulation, and pacemaker interference.

#### 5.3.1. Monitor (Display) Interference

##### 5.3.1.1. Continuous Interference

1. Check the Power Cord connections for the ASU.
2. Check all other electrical equipment in the operation room for defective grounds.
3. If the electrical equipment is grounded to different objects, rather than a common ground, voltage differences can appear between the two grounded objects. The monitor may respond to these voltages. Some types of input amplifiers can be balanced to achieve optimum common mode rejection and may possibly correct the problem.

##### 5.3.1.2. Interference Only When ASU is Activated

1. Check all connections to the ASU, and active accessory to look for possible metal-to-metal sparking.

- If interference continues when the ASU is activated and while the electrode is not in contact with the patient, the monitor is responding to radio frequencies. Some manufacturers offer RF choke filters for use in the monitor leads. These filters reduce interference while a generator is activated. RF filters minimize the potential for an electrosurgical burn at the site of the monitor electrode.
- Check that the ground wires in the operating room are electrically consistent. All ground wires should go to the same grounded metal with wires that are as short as possible.
- If the above steps do not remedy the situation, qualified service personnel should check the ASU.

### 5.3.2. Neuromuscular Stimulation

- Stop the surgery.
- Check all connections to the ASU and active electrodes to look for a possible metal-to-metal spark.
- If no problems are found, the ASU should be checked by qualified service personnel for abnormal 50/60 Hz AC leakage current.

### 5.3.3. Pacemaker Interference

- Check all connections.
- Always monitor pacemaker patients during surgery.
- Always keep a defibrillator available during electrosurgery on patients with pacemakers.
- Consult the pacemaker manufacturer for specific recommendations.

## 6. SYMBOLS USED

	Alternating Current		Volume Control		Caution		Non-Sterile
	Caution: Electrical Shock Hazard		Manufacturer		Type CF Defibrillation-Proof Applied Part		Consult Instructions for Use
	READY		RF ON		Transmurality		Equipotential Terminal
	Footswitch		Catalogue Number		Serial Number		Waste Electrical and Electronic Equipment (WEEE)
	Fuse Information		Model Number		Non-Ionizing Radiation		Authorized representative in the Brazilian Community
	Country and Date of Manufacture		Does not contain phthalates		Medical Device		Not Made with natural Rubber Latex
	Input Voltage Selector Switch		Data Port		Unique Device Identifier		Keep Upright
	Power ON		Power OFF	<b>Rx ONLY</b>	Prescriptive Use Only		
	Service Access		Dangerous Voltage			National Institute of Metrology Standardization and Industrial Quality	
		Humidity Limit (Storage and Transit)					
		Temperature Limit (Storage and Transit)					

## 7. TECHNICAL SPECIFICATIONS AND SAFETY INSPECTION

### 7.1. RF Output

- Frequency: 460 kHz  $\pm$ 5%, Quasi-sinusoidal
- ASU Maximum Power Output: 32.5 W at 100 $\Omega$
- HF Power and Voltage Output:

Device Code	Maximum Output Power	Maximum Output Voltage	Handpiece Type
A	28.5 W at 114 $\Omega$	57.0 Vrms	Isolator <sup>®</sup> Synergy™ clamp
B	15.0 W from 20 $\Omega$ to 400 $\Omega$	77.5 Vrms	Isolator™ Transpolar™ pen
C	20.0 W from 31 $\Omega$ to 300 $\Omega$	77.5 Vrms	Isolator™ Transpolar™ pen Isolator™ linear pen
D	25.6 W at 127 $\Omega$	57.0 Vrms	Isolator™ clamp
E	22.8 W at 143 $\Omega$	57.0 Vrms	Isolator™ clamp
F	28.5 W at 114 $\Omega$	57.0 Vrms	Isolator™ clamp
G	28.5 W at 114 $\Omega$	57.0 Vrms	Isolator™ clamp
H	28.5 W at 114 $\Omega$	57.0 Vrms	Isolator™ clamp
J	12.0 W from 20 $\Omega$ to 500 $\Omega$	77.5 Vrms	Isolator™ Transpolar™ pen
K	25.0 W from 39 $\Omega$ to 240 $\Omega$	77.5 Vrms	Isolator™ Transpolar™ pen
L	30.0 W from 47 $\Omega$ to 200 $\Omega$	77.5 Vrms	Isolator™ Transpolar™ pen

### 7.2. Mechanical Specifications and Software Version

#### Size:

**ASU:** 17.5" (44.5 cm) D x 13.75" (35 cm) W x 6" (15 cm) H maximum.

#### Weight:

**ASU:** 15 lbs. (6.8 kg) absolute maximum; 20 lbs. (9 kg) including cables

#### Software Version

V3.11B

### 7.3. Environmental Specifications

Operating Conditions	
Temperature	10°C to 40°C, 50°F to 104°F
Humidity	15% RH to 90% RH, non-condensing
Atmospheric pressure	80 kPa to 106 kPa, 11.603 psi to 15.374 psi
Storage Conditions	
Temperature	-29°C to 60°C, -20°F to 140°F
Humidity	30% RH to 85% RH, non-condensing
Transit Conditions	
Temperature	-29°C to 60°C, -20°F to 140°F
Humidity	30% RH to 85% RH, non-condensing

### 7.4. Electrical Specifications

- 100-120V ~ 50/60 Hz

### 7.5. Fuse Specifications

- Littelfuse 1.25A 250V, T-lag, 5 x 20 mm, UL Recognized, IEC Approved, RoHS

**⚠ Caution:** Do not alter the fuse from factory settings.

### 7.6. Footswitch Specifications

- Moisture protection rating: **IPX8**

## 7.7. Power and Voltage Output Restrictions

The maximum power output of 28.5 W for the Isolator® Synergy™ clamp is available at 114Ω load for devices operating under device code “A, F, G, and H”. Lower maximum output powers are available depending on the system operating mode. See Section 7.1.

The maximum power output of 15.0 W for the Isolator™ Transpolar™ pen is available between 40Ω to 400Ω load for device operating under device code “B”. Lower maximum output powers are available depending on the system operating mode. See Section 7.1.

The maximum power output of 20.0 W for the Isolator™ linear pen is available between 31Ω to 300Ω load for device operating under device code “C”. Lower maximum output powers are available depending on the system operating mode. See Section 7.1.

At other load impedances, the ASU will reduce the available power to comply with the specified voltage and current limits. See Figure 11 and Figure 12.

The ASU is capable of producing a maximum output power of 32.5 Watts under a 100 Ohm load although no current AtriCure® Bipolar Handpiece uses power above 30 Watts.

The maximum output voltage depends on the device code and can be either 57 Vrms or 77.5 Vrms. See Section 7.1.

## 7.8. Equipment Type / Classification

- Class 1 Equipment


Figure 11 – Power vs. Load (clamp algorithm)


Figure 12 – Power vs. Load (pen algorithm)

## 8. PREVENTIVE MAINTENANCE AND CLEANING OF ASU

### 8.1. Preventive Maintenance

Perform annual preventative maintenance procedures to ensure all ASU components are functioning as defined within this manual. Pay particular attention to operational and safety features, including but not limited to:

- Electrical power cords for fraying, damage, and proper grounding
- AC power switch
- Indicator damage (Power On, Fault, Ready, RF ON, Transmurality)
- LCD display damage or loss of graphic information
- Temperature display damage or loss of information
- Handpiece connector damage, cracking or inability to insert and latch Handpiece plug
- Carrying handle damage, inability to latch or rotate
- Rubber feet damage, cracking or inability for the ASU to remain stable on a flat surface.
- Footswitch cord fraying or damage
- Footswitch connector damage cracking or inability to insert and latch footswitch plug
- Footswitch pedal damage check activation by pressing and releasing the pedal

Other medical equipment that may be used simultaneously with the ASU should also be inspected for damage. Specifically, check for insulation damage of monitoring electrode cables and endoscopically used accessories.

Visually inspect the footswitch for fluids or other infectious hazards. Clean as necessary using the instructions in Section 8.2.

The ASU does not have any serviceable parts. For servicing issues, contact AtriCure, Inc. at:

AtriCure, Inc.  
7555 Innovation Way  
Mason, Ohio 45040  
USA  
+1 866 349 2342  
+1 513 755 4100

### 8.2. Cleaning and Disinfecting

**NOTE: Do not spray or pour liquids directly on the unit.**

**NOTE: The unit and/or components cannot be sterilized.**

#### WARNING

70% Isopropyl Alcohol (IPA) alcohol wipes should be used to clean the exterior of the ASU/ASB. Do not submerge or allow fluids to enter the chassis. Do not spill fluid over the ASU/ASB system and ensure IPA is completely dry before operating the units to avoid damage to the equipment or injury to the patient. If the fluid is spilled on the ASU/ASB, return it to the hospital's Biomedical Engineering Department for evaluation.

 **CAUTION:** Avoid caustic or abrasive cleaners

 **CAUTION:** If the tamper proof seal is broken, the integrity of the ASU/ASB cannot be ensured. Contact AtriCure representative.

#### Guidelines

The following guidelines are recommended for cleaning the unit. It is the user's responsibility to qualify any deviations from these processing methods.

1. Disconnect the unit or cart from the outlet before cleaning.
2. If the unit and/or components are contaminated with blood or other body fluids, they shall be cleaned before the contamination can dry (within two hours of contamination).
3. The outer surfaces of the unit and/or components shall be cleaned with 70% Isopropyl alcohol (IPA) wipes for a minimum of two minutes. Do not allow fluids to enter the chassis.
4. Pay attention to all areas where fluids or soil may gather, such as under/ around the handles or any tight crevices/ grooves.
5. Dry the unit and/or components with a dry, white lint-free cloth.
6. Conduct a final confirmation of the cleaning process by visually inspecting the white cloth for remaining soil.
7. If soil remains on the white cloth, repeat steps 3 through 6.
8. After cleaning is complete, turn the unit on to perform Power On Self-Test (POST). If any errors are received, contact AtriCure to begin return process.

**⚠️ WARNING ⚠️**

If necessary, clean the RF Interface Cable, Footswitch Interface Cable, and the PSS Cables with 70% IPA wipes to ensure it does not cause infection. Inappropriate handling of the Interface or PSS cables, including sterilization and immersing the electrical connectors, can result in degradation of system performance, including the inability to initiate or complete an ablation therapy.

**⚠️ CAUTION:** Ensure all Interface cables and Footswitch are adequately cleaned/disinfected.”

**⚠️ CAUTION:** Ensure IPA is completely dry before operating the cables.

## 9. DISPOSAL

Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

## 10. ACCESSORIES

### 10.1. ASB3, Switch Matrix Accessory

The Switch Matrix provides a means of connecting multiple handpieces to the ASU and a means of selecting the input to the handpiece electrodes. The input is selected with the Switch Matrix knob. A cable with is provided to connect the Switch Matrix to the ASU.

**⚠️ WARNING ⚠️**

**Do not connect the ASB3 auxiliary device cable to supply mains (line voltage) operated equipment without evidence that the safety certification of the Accessory has been performed in accordance to the appropriate EN60601-1 and/or EN60601-1-2 harmonized national standard. Supply mains operated equipment may introduce dangerous leakage currents into the heart.**

An auxiliary device (other than those listed in paragraph 10.2.2) may have an adverse effect on nearby radio or TV or medical equipment. There may also be cases when nearby electrical appliances adversely influence the auxiliary device, causing data errors or malfunction.

Auxiliary devices compatible for use with the Switch Matrix include:

- Any AtriCure Isolator™ Handpiece
- Any AtriCure Transpolar™ Pen

**⚠️ WARNING ⚠️**

Read auxiliary device manual and observe warnings.


Any AtriCure handpiece or pen device may be connected to the Switch Matrix. AtriCure Devices will be functional when the device is connected to the correct receptacle and the Switch Matrix switch knob is turned to indicate the device for use.

**⚠️ Caution:** Ensure the ASB Selector Switch is positioned on the correct Handpiece electrode input.

**⚠️ Caution:** Ensure Handpieces and cables are fully and properly installed into the correct receptacle.

Settings and procedures for the auxiliary device are determined according to the instructions for use provided with the auxiliary device.

ASB3, Switch Matrix Unit set up is shown below with the ASU.

## 11. COMPONENTS AND CABLES

- ASU/ASB Power Cord
- ASU/ASB RF Interface Cable
- ASU Footswitch
- ASU/ASB Footswitch Interface Cable

# WARRANTIES

## Limitation on Liability

This warranty and the rights and obligations hereunder shall be construed under and governed by the laws of the State of Ohio, U.S.A.

AtriCure, Inc. warrants this product to be free from defects in material and workmanship under normal use and preventive maintenance for the respective warranty period shown below. AtriCure's obligation under this warranty is limited to the repair or replacement, at its option, of any product, or part thereof, which has been returned to AtriCure, Inc. or its Distributor within the applicable time period shown below and which examination disclosed, to AtriCure's satisfaction, to be defective. This warranty does not apply to any product, or part thereof, that has been: (1) adversely affected due to use with devices manufactured or distributed by parties not authorized by AtriCure, Inc. (2) repaired or altered outside AtriCure's factory in a way so as to, in AtriCure's judgment, affect its stability or reliability, (3) subjected to improper use, negligence or accident, or (4) used other than in accordance with the design and use parameters, instructions and guidelines for the product or with functional, operational or environmental standards for similar products generally accepted in the industry. **AtriCure has no control over the operation, inspection, maintenance or use of its products after sale, lease or transfer, and has no control of the selection of Customer's patients.**

AtriCure's products are warranted for the following periods after shipment to the original purchaser:

ATRICURE ABLATION AND SENSING UNIT.....	ONE (1) YEAR
ATRICURE SWITCH MATRIX.....	ONE (1) YEAR
RF AND FOOTSWITCH INTERFACE CABLES.....	ONE (1) YEAR
ATRICURE FOOTSWITCH.....	ONE (1) YEAR
GROUNDED ELECTRICAL CORD.....	ONE (1) YEAR

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ATRICURE, INC. AND IS A PURCHASER'S EXCLUSIVE REMEDY. IN NO EVENT SHALL ATRICURE, INC. BE LIABLE FOR SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES INCLUDING, WITHOUT LIMITATION, DAMAGES RESULTING FROM LOSS OF USE, PROFITS, BUSINESS OR GOODWILL.

AtriCure, Inc. neither assumes nor authorizes any other person to assume for it any other liability in connection with the sale or use of any of AtriCure Inc. products. There are no warranties that extend beyond the terms presented unless an extended warranty is purchased before the original warranty expires. **No agent, employee or representative of AtriCure has any authority to change any of the foregoing or assume or bind AtriCure to any additional liability or responsibility.** AtriCure, Inc. reserves the right to make changes to products built and/or sold by them at any time without incurring any obligation to make the same or similar changes on products previously built and/or sold by them.


### DISCLAIMER

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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