AtriClip® Selection Guide

Instructions for Use

CGG100

MD

CAUTION: Federal law (US) restricts this device to sale by or on the order of a physician.

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AtriClip® Selection Guide

INDICATION FOR USE

The AtriClip Selection Guide (Guide) is intended to be used as an accessory to the AtriClip LAA Exclusion System with preloaded Clip (Clip).

CONTRAINDICATIONS

None known

DESCRIPTION

The Guide is an accessory that works in conjunction with the AtriClip LAA Exclusion System. The AtriClip LAA Exclusion System is comprised of the AtriClip LAA Exclusion Device (Clip), the AtriClip Applier, and the Guide.

The Guide is a single patient use surgical instrument designed to assist in the selection of the appropriate Clip. The Guide is malleable and may be placed directly adjacent to the Left Atrial Appendage (LAA) for sizing. The markings on the Guide facilitate the evaluation of the structure and selection of the appropriate Clip size. The markings on the Guide are 4 mm +/- 0.5 mm (0.16 in +/- 0.02 in) shorter than the Clip nominal dimensions to approximate the appendage compression when the Clip is applied. Tissue thickness, appendage geometry, and other factors can affect sizing decisions. Therefore, it is the physician's judgment to select the appropriate size.

The Guide (CGG100) is a sterile accessory composed of aluminum and cured polyurethane ink. It is not made with natural rubber latex and does not contain phthalates.

ENVIRONMENTAL SPECIFICATIONS

Storage	Transit
Temperature: -29°C/-20°F to 60°C/140°F	Temperature: -29°C/-20°F to 60°C/140°F
Relative Humidity: 15% to 85%	Relative Humidity: 30% to 85%
Atmospheric Pressure: N/A	Atmospheric Pressure: N/A

PACKAGE CONTENTS

- One (1) Guide

The Guide is supplied STERILE and NON-PYROGENIC in an unopened, undamaged package. For single use only. Do not re-sterilize. Do not re-use.

ATRICLIP SELECTION GUIDE

NOMENCLATURE (SEE FIGURE 1)

- Handle
- Base
- Indication Mark

⚠ WARNINGS ⚠

Do not apply excessive force when using the Guide. Using excessive force may cause tissue damage.

Read all instructions for the Guide before use and use the device only as intended. Use of the Guide should be limited to properly trained and qualified medical personnel. Improper use of this device may lead to device malfunction, failure to provide intended therapy, and/or serious injury.

Do not bend the Guide in the area of the indication marks. This may lead to incorrectly determining corresponding AtriClip size. Failure to correctly size the Clip may result in: tissue trauma, dehiscence, tissue tearing, displacement, and/or lack of desired hemostasis.

DO NOT RESTERILIZE. The Guide is provided STERILE and is intended for SINGLE use only. Re-sterilization may cause injury to the patient.

Use caution when using the Guide to determine the corresponding AtriClip size. Failure to correctly size the Clip may result in: tissue trauma, dehiscence, tissue tearing, displacement, and/or lack of desired hemostasis.

INSTRUCTIONS FOR USE

- Using sterile technique, remove the Guide from its packaging.

⚠ WARNING ⚠

If the sterile package is damaged and/or the sterile barrier is breached, discard the device and DO NOT USE to avoid the risk of patient infection.

CAUTION: Do not drop the device as this may induce damage to the device. If the device is dropped, do not use. Replace with a new device.

- Align the Guide adjacent with the LAA (See Figure 2).

NOTE: Evaluation should occur at the base of the appendage.

- Evaluate the LAA and select the appropriate Clip size.

CAUTION: The Guide is to be used only to assist with selecting the appropriate Clip.

Labeled Clip Size	LAA Size Range
35 mm	29 - 35 mm (1.14 - 1.38 in)
40 mm	34 - 40 mm (1.34 - 1.57 in)
45 mm	39 - 45 mm (1.54 - 1.77 in)
50 mm	44 - 50 mm (1.73 - 1.97 in)

DISPOSAL INFORMATION

After use, this device should be treated as medical waste and disposed of following hospital protocol.

SERIOUS INCIDENT

Any serious incident that has occurred in relation to this device should be reported to AtriCure.

RETURN OF USED PRODUCT

If for any reason this product must be returned to AtriCure, Inc., a return goods authorization (RGA) number is required from AtriCure, Inc., prior to shipping.

If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from AtriCure, Inc.

DISCLAIMER STATEMENTS

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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AtriCure Inc.