AtriCure®

EPI-Sense® ST RF Cable (“RF Cable”) Instructions for Use

CSK-2060

MD

INSTRUCTIONS FOR USE

EPI-Sense® ST RF Cable

DEVICE DESCRIPTION:

Sterile cable for connection of compatible EPI-Sense® ST device (EPI-Sense ST Coagulation Device) and compatible AtriCure RF Generator (CSK-320 or MAG) or Sensing Cable. STERILE in unopened, undamaged package. For single use only. Do not re-sterilize. Do Not Re-Use.

CONTENTS:

EPI-Sense ST RF Cable contains:

- One (1) CSK-2060 EPI-Sense ST RF cable in sterile packaging

The EPI-Sense ST RF Cable is not made with natural rubber latex and does not contain PVC or phthalates.

INDICATIONS FOR USE:

The RF cable is designed to transmit RF energy from the compatible AtriCure RF Generator to the attached EPI-Sense family coagulation device to achieve its intended purpose.

INTENDED PURPOSE:

The RF cable is designed to transmit RF energy from the compatible AtriCure RF Generator to the attached EPI-Sense family coagulation device to achieve its intended purpose.

INTENDED USER AND PATIENT POPULATION:

The RF cable is a medical device for use by licensed medical doctors who perform cardiac and/or thoracic procedures using the EPI-Sense family of coagulation devices for treatment of patients with longstanding persistent atrial fibrillation.

INTENDED CLINICAL BENEFIT:

To achieve the clinical benefits of the attached EPI-Sense family of coagulation devices.

⚠ WARNING ⚠
For information concerning all warnings, precautions, and troubleshooting, see compatible AtriCure RF Generator (CSK-320 or MAG) Operators Manual and attached EPI-Sense family coagulation device IFU. Failure to follow the instructions contained in the RF generator manual or coagulation device IFU may lead to an inability to complete the procedure.
For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
Inspect all devices and packaging prior to use. If any breach of the packaging is found the sterility of the product cannot be ensured which poses a risk of patient injury. Do not use product if breach is found.
Do not use if cable is damaged as the procedure may not be able to be completed. Replace cable.
Do not use cable if yellow LED on EPI-Sense family coagulation device is illuminated. This indicates low battery and RF ablation energy will not be delivered. Replace cable. The RF cable has a limited functional life; if greater than 4 hours of use are required, it is recommended to use a second cable.
The Coagulation Device, compatible AtriCure RF Generator, Cables and Accessories have been tested as a system and comply with the limits for medical devices to IEC 60601-1-2. Use of another manufacturer’s accessories may cause damage to the equipment or injury to the patient.
Replace cable if additional ablations cannot be performed after 30 ablations.

INSTRUCTIONS FOR USE:

- Inspect package and RF cable for damage. If either RF cable or package is damaged, discard without using.
- Remove RF cable from sterile package.

⚠ WARNING ⚠
Arrow indicators and color coded ends on connectors should be used to aid in correct alignment to avoid damage to cable ends.

⚠ CAUTION: Cables to surgical electrodes should be positioned to prevent contact with patient or other leads.

- Connect BLUE strain relief end of cable to BLUE strain relief end of the EPI-Sense family coagulation device.
- Connect BLACK end of cable to:
 - BLACK port on RF Generator (CSK-320 or MAG) or
 - BLACK port on the sensing cable
- Follow instructions for use of generator, coagulation device, and sensing cable.
- Do not disconnect cable, grasp by connector cover and pull back. Do NOT pull on cable.

⚠ WARNING ⚠
RF Cable contains a Lithium disposable battery. Do not recharge, disassemble, heat above 60°C, incinerate, or expose the battery directly to water to avoid injury to user or damage to device.
Ensure device is disposed of following local governing ordinances and recycling plans to prevent biohazard exposure.

- After use, this device should be treated as medical waste and disposed of following hospital procedure.
 - To remove battery, unscrew the cover, using a #1 Phillips screwdriver, to access the battery.

	Temperature	Humidity	Atmospheric Pressure
Operational	10° to 40°C (50° to 104°F)	30%-75% relative humidity non-condensing	700 to 1060 millibar (10 to 15 psi)
Transit	-20° to 55°C (-4° to 131°F)	30-85% relative humidity	N/A
Storage	-20° to 55°C (-4° to 131°F)	30-85% relative humidity	N/A

DISPOSAL

After use this device should be treated as medical waste and disposed of following hospital protocol.

RETURN OF USED PRODUCT

If for any reason this product must be returned to AtriCure, Inc., a return goods authorization (RGA) number is required from AtriCure, Inc., prior to shipping. If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from AtriCure, Inc.

NOMENCLATURE:

- RF = Radiofrequency
- IFU = Instructions For Use
- LT = Label

LIMITED WARRANTY:

AtriCure warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. AtriCure’s sole obligation under this warranty is limited to the repair or replacement of this instrument. AtriCure neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond AtriCure’s control directly affect the instrument and the result obtained from its use. AtriCure assumes no liability with respect to instruments deliberately mis-used or those reused, reprocessed or re-sterilized and makes no warranties expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such mis-used or reused instruments. AtriCure shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the deliberate mis-use or re-use of this instrument.

DISCLAIMER:

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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